

Multi-type Optional

3 Combinations

- Logical bending + 9 o'clock instrument exit
- Anti-Logical bending + 3 o'clock instrument exit
- Anti-Logical bending + 9 o'clock instrument exit



Ordering information

Product	Item code
Single-use ureterorenoscope-HU30M	
EU Standard (Anti-logic Deflection) (L)	B-8000-043
US Standard (Logic Deflection) (L)	B-8000-053
EU Standard (Anti-logic Deflection) (R)	B-8000-052
Medical Image Processor	
HUV-02 with its accessories	2201B-PA001
MS-8 with its accessories	2288B-PA001

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

HugeMed

f-URS enters the era of

6.3_{FR}

Single-use
Ureterorenoscope

HU30M



HugeMed

Shenzhen HugeMed Medical Technical Development Co., Ltd.
401, 501, Building4, FutongHaizhi TechnologyPark, No.17, BulanRoad, XialilangCommunity, NanwanStreet,
LonggangDistrict, Shenzhen, Guangdong, China, 518112


+86-755-22275833 +86-755-22275833 info@hugemed.net

For more information
please visit our website:
www.hugemed.net



P/N: ENG-EU-CATA-fURS-suHU-HU30M-202404V1.0

HugeMed HU series f-URS EVOLUTION







- The cross-section of HugeMed™ f-URS
- The cross-section of the working channel
- The cross-section of the expulsion space
- The cross-section of the UAS
- TESIR: Theoretical expulsion space increasing rate

**Note: The comparison of the expulsion space indicates:
HU30 < (39%) HU30S < (27%) HU30M.**

What clinical challenges
are addressed by

SLIMMER





Enhanced Stone Clearance Capacity

- Shorten RIRS lithotripsy operation time
- Reduce incidence of post-RIRS Stone Street
- More complete stone-free rate of single-session RIRS

Decrease Stenting


Procedure in pre-operation and post-operation

- More primary ureteroscopy
- Reduce patient length of stay
- An alternative to emergency stenting



Offer More Space to Intra-Renal Temperature and Pressure Management

- More outflow space for taking away laser heat
- More outflow space for releasing irrigation pressure



What else has been
Upgraded



Connection
Gold Fingers
PCB plug

- Smooth insertion
- Stable electrical conductivity




Flexibility

120° left and right
Insertion tube rotation

- Added mobility of exploration
- More accurate to aim at targets
- Enhanced team collaboration


Compatible With 2 Types of Processors



HUV-02

Medical Imaging Processor

- Touch-screen control
- Standard edge connector
- Flat design for easy deployment
- Full HD output 1920*1080 resolution



MS-8

Medical Imaging Processor

All-in-One

Display | Processor | Recorder
15.6 inch touch control

- A portable imaging solution for ureteroscopy clinical practice.